

Viva Zym

Concentrated enzyme digester

Viva Zym is a highly active enzymatic concentrate for removal of all blood and albuminous stains and soiling. Best results are obtained at 100 to 140 °F and a pH-value below 10.

Physico-chemical data:

Consistency:	homogeneous liquid	
Total alkalinity:	n.d.	(1 % solution)
Active alkalinity:	n.d.	(1 % solution)
pH-value:	~ 8.0	(1 % solution)
Density:	~ 1.08 oz/fl.oz	(68 °F)

Recommended Dosage:

Per kg dry laundry at a bath ratio of 1 : 4

Hardness Range	Water Hardness mmol CaCO ₃ / l	1-Bath Procedure	2-Bath Procedure	
			1 st Bath	2 nd Bath
1	< 1,5 soft	0.1 – 0.8 fl.oz	0.1 – 0.8 fl.oz	
2	1,5 – 2,5 medium	0.1 – 0.8 fl.oz	0.1 – 0.8 fl.oz	
3	> 2,5 hard	0.1 – 0.8 fl.oz	0.1 – 0.8 fl.oz	

Ecology:

Viva Zym is exempt from chlorine, APEO, EDTA, NTA, and is easily biodegradable.

Storage conditions: To be stored in a cool and dry place, protected from frost. Protect from heat and direct sunlight.

Maximum storage : 12 months

HEALTH AND SAFETY RECOMMENDATIONS:

PRECAUTIONS: Keep out of the reach of children. When using, do not eat or drink. Avoid contact with skin and eyes.

FIRST AID: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, take off immediately all contaminated clothing and wash immediately with plenty of water. In case of accident or if you feel ill, seek medical advice immediately. If swallowed, seek medical advice immediately and show this container or label.

ONLY TO BE USED IN PROFESSIONAL LAUNDRY SYSTEMS!

FOR INDUSTRIAL USE ONLY!

The information provided herein is given on an advisory basis which does not imply any liability on our part. For further details see our relevant safety data sheet prior to use !